OPTIFIX®

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Operating manual

Mode d'emploi

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Hersteller/Manufacturer:
Poulten & Graf GmbH
Am Bildacker 3 - 7
97877 Wertheim, Germany

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

**HF
SAFETY S**

Preface

We thank you for the purchase of an OPTIFIX® bottle top dispenser. You made a good choice. With more than 20 years of experience in making dispensers we guarantee a high standard of product quality and safety. To guarantee quality and safety of the dispenser for the long term, proper operation and maintenance is necessary. Therefore please read this operating manual carefully before using the dispenser and take note of the warning instructions as well as the corresponding regulations. The pack contains: 1 dispenser (valve block unit with finger guard), 1 discharge tube set (model SAFETY, SAFETY S and HF 1 additional stopper cap), 1 intake tube, 2 adapters (model BASIC and SOLVENT: made of PP; model SAFETY, SAFETY S and HF: made of PTFE), 1 operating manual. Please check this before using the dispenser.

Description of Cross-section Drawing

1. Valve block made of solid PTFE.
2. FIX adaptor with internal thread GL 32 made of PP (2-10 ml) from 30 -500 ml + HF: PTFE)
3. Ventilation channel with BASIC, SOLVENT and SAFETY.
- 3a. Ventilation channel and ventilation valve with SAFETY S and HF.
4. Intake valve with sapphire ball.
5. Discharge valve with sapphire ball.
6. Discharge tube set (consists of: PTFE-discharge tube; support made of PVDF, 100 ml made of PP, 200 - 500 ml: made of stainless steel; PTFE-thread nipple; PVDF-screw cap; PP-stopper cap).
7. Finger guard with integrated scale made of PP.
8. Volume fine adjustment screw made of PP.
9. Volume setting system made of PP.
10. Stopper cap made of PP.
11. Glass cylinder, plastic coated with stopper ring.
12. Glass piston, solid PTFE coated, patented; additional seal at the bottom of the piston with model SAFETY, SAFETY S and HF. Model SOLVENT with precision glass piston.

Operating Manual Bottle Top Dispenser OPTIFIX®

I. Safety Instructions

1. Before using the dispenser please read the operating manual carefully and follow the safety instructions. **This operating manual cannot reasonably mention all safety risks that might be associated with the use of the dispenser. It is therefore in the responsibility of the user to consult where necessary and establish appropriate health and safety practices to determine the applicability of regulatory limitations prior to use.**
2. Observe general notes on hazards and safety regulations, e. g. wear protective clothing, gloves and glasses.
3. Observe any hazard symbols on reagent bottles.
4. Use only appropriate containers to dispense liquids. If using small bottles, please use a stand.
5. Always make sure that the discharge tube never points at the user or any other person.
6. When transporting the dispenser please ensure that the piston is in zero-position and no liquid is in the cylinder.

The stopper cap must close the discharge tube and the volume setting system must be set at zero-position! Never carry the dispenser by the finger guard, only by the PTFE-valve block!

7. Before use, please check correct assembly of the dispenser with the bottle, that the bottle is properly fitted and closed to full lock. This will be possible using standard glass laboratory bottles with a standard thread (e.g. Schott bottles). By use of short neck bottles, bottles or containers with a short thread, or other bottles where full tightening of the bottle top dispenser onto the bottle cannot be achieved, the dispenser should be screwed onto the bottle using an adaptor to ensure secure operation. Most adaptors are readily available from our catalogue or can be supplied according to special requirements on request. **The stopper cap of the discharge tube may never come into the dosing jet, in order to avoid dangerous squirting.**

8. To guarantee trouble-free working please clean and use the dispenser according to the operating manual.

9. Repairs if necessary must be carried out only by the manufacturer or authorized service stations using original parts.

10. For safety reasons we recommend use of adaptors made of PTFE which are available as an option.

The models SAFETY, SAFETY S and HF must be used with adaptors made of PTFE.

II. The World of OPTIFIX®

OPTIFIX® BASIC (Art. No. 101.080)

Bottle top dispenser for aqueous liquids, acids, bases and crystallizing liquids. Piston made of glass, solid PTFE coated, patented. Safety dispensing cylinder, plastic coated. Valve block made of pure PTFE. Valve system without metal springs. FIX adaptor.

OPTIFIX® SOLVENT (Art. No. 101.081)

Bottle top dispenser for solvents. Piston made of precision glass. Safety dispensing cylinder, plastic coated. Valve block made of pure PTFE. Valve system without metal springs. FIX adaptor.

FORTUNA UNIVERSAL® OPTIFIX® HF (Art. No. 101.106)

Bottle top dispenser for extremely hazardous chemicals, especially hydrofluoric acid. Piston made of glass, solid PTFE coated (patented), with seal. Dispensing cylinder and valve block made of pure PTFE. Valve system without metal springs. Ventilation with ball valve. FIX adaptor.

OPTIFIX® SAFETY (Art. No. 101.107)

Bottle top dispenser for strong acids or hazardous chemicals. Piston made of glass, solid PTFE coated (patented), with seal. Safety dispensing cylinder, plastic coated. Valve block made of pure PTFE. Valve system without metal springs. FIX adaptor.

OPTIFIX® SAFETY S (Art. No. 101.107S)

Bottle top dispenser for strong or fuming acids or highly odorous chemicals. Piston made of glass, solid PTFE coated (patented), with seal. Safety dispensing cylinder, plastic coated. Valve block made of pure PTFE. Valve system without metal springs. Ventilation with ball valve. FIX adaptor.

Never use the OPTIFIX® Dispensers with:

- liquids attacking Al_2O_3 and PTFE
- liquids attacking borosilicate glass (except with the OPTIFIX® HF)

Warning:

The BASIC and SOLVENT model come with adaptors made of polypropylene to fit on most popular containers. However, when dosing extremely hazardous chemicals, for safety reasons it is recommended to use adaptors made of PTFE for BASIC and SOLVENT model which are available as an option.

III. Modern Technique

Finger guard made of PP. All descriptions on the finger guard and the graduation on the scale are lasered to guarantee permanence.

Volume setting system Easy-to-use volume setting system with fine adjustment for calibration. Suitable for use by left- and right-handers.

Dispensing cylinder Glass cylinder, vacuum-calibrated, plastic coated, with stopper ring. Black graduation. Model HF with dispensing cylinder made of pure PTFE.

Valve block made of pure solid PTFE, with FIX Adaptor.

Valve seats made of glass, precision ground for precision sapphire valve balls. No metal springs. Works by means of gravity. Model HF with valve seats made of sapphire.

Dispensing piston Model BASIC, SAFETY, SAFETY S and HF glass piston, coated with solid PTFE, patented. Seal with model SAFETY, SAFETY S and HF. Model SOLVENT with glass piston.

Discharge tube set consists of: discharge tube, support, screw cap, thread adaptor and stopper cap.

IV. Handling


1. Assembling / Filling

Before using the dispenser please check the dispenser for any damage or glass breakage.

Push intake tube over the hose connection of the PTFE valve block.

Screw dispenser direct onto containers with thread GL 32 or use adaptors to fix on any other threads or ground joint necks. Check the correct fitting. The dispenser can be turned on the flask and fixed in any direction by using the built-in FIX adaptor.


Screw the discharge tube set onto the thread adaptor of the PTFE valve block. Remove stopper cap from the discharge tube.

Warning:


If the stopper cap is not removed when dispensing, there may be an overpressure in the channel system and uncontrolled dispensing may be the result.

To fill the dispenser the first time lift and lower piston a short distance several times until no air bubbles are seen in the cylinder. Model HF: Because of the PTFE-cylinder there are no graduations on the cylinder. (Note: - bleeding air from the system with return of all liquid direct to the reservoir is possible with the 'LOOP' system)

2. Volume setting

Loosen the nut of the volume setting system, set the desired volume. The upper edge of the fine adjustment screw shows the volume set. Fix the screw.

For calibrating purposes turn the screw of the fine adjustment. Turn downwards to increase the volume, turn upwards to decrease the volume.

3. Dispensing procedure

- a) For safety make sure that the discharge tube does not point at the user or any other person. By using the built-in FIX adapter it is possible to fit the dispenser as required. Remove stopper cap by holding the discharge tube tight with other hand.
- b) Set the desired volume with the volume setting system.
- c) Lift the dispensing piston up with the finger guard. Do not touch the piston itself to avoid any unlikely contact with chemicals. The piston stops at the volume set. The graduation on the glass cylinder is an additional optical control of the volume set or air-bubbles. There is no graduation on the PTFE-cylinder of the model HF.

d) The models BASIC and SOLVENT work by gravity. To dispense liquid let the piston go; it will fall by its own weight. Do not press by hand!

The models SAFETY, SAFETY S and HF are equipped with a special piston seal. To dispense liquid push the piston down to the zero-position using the finger guard. **The stopper cap of the discharge tube may never come into the dosing jet, in order to avoid dangerous squirting.**

e) Avoid uneven piston movements or jerking of the stopper ring.

f) Replace stopper cap after dispensing by holding the discharge tube tight with other hand.

Warning:

The piston of the models BASIC and SOLVENT run down to the zero-position automatically under their own weight. With the models SAFETY, SAFETY S and HF with special seal the piston must be pushed down manually. Therefore chemicals may remain in the dispensing cylinder. To avoid unintentional dispensing, stopper cap must be placed on the discharge tube after every use. When dispensing observe safety precautions!

V. Operating Temperature

The operating temperature for the models BASIC and SOLVENT is from +15° C (+59° F) to +50° C (+122° F), with the models SAFETY, SAFETY S and HF from +20° C (+68° F) to +30° C (+86° F).

VI. Sterilization

The models BASIC and SOLVENT can be sterilized at +121°C (+248°F)/1 bar. The dispensing piston of the SOLVENT can remain in the cylinder during the sterilization, the BASIC piston has to be sterilized separately.

Before sterilisation loosen screw of the volume setting system. **Attention:** The BASIC sizes 200, 300 and 500 ml are not prepared for sterilization but they are available in a modified version suitable for sterilization. Please ask your dealer or the manufacturer.

The models SAFETY, SAFETY S and HF cannot be sterilized. Please clean with appropriate cleaning solutions or ultrasound.

Sterilisation of OPTIFIX® dispensers will put excessive stress on the dispensers both thermally and mechanically. As a consequence sterilisation of OPTIFIX dispensers will invalidate the warranty offered on the dispensers.

VII. Transport and Storage

1. When moving the dispenser never carry by the finger guard but always by the PTFE- valve block to avoid unintentional dispensing of liquid! Replace the stopper cap before transporting and fix the volume setting system at the zero-position.

2. Please transport and store the dispenser only in the original box to avoid damage. Storage temperature: -10° C (+14° F) to +50° C (+122° F).

VIII. Maintenance and Cleaning

1. To guarantee trouble-free use for many years it is necessary to clean the dispenser regularly. **If there is any sign of a potential malfunction (e.g. piston is difficult to move) never use force. Immediate cleaning should be carried out.**
2. Do not leave the dispenser empty without cleaning to avoid clogging the valve balls! If this does happen soak the dispenser for some hours in cold water or clean it using ultrasound. Once the piston can be moved, dismantle and remove any residual impurities by cleaning with a normal solution.
3. The dispenser can be left on the bottle providing the channel system is filled.
4. When changing the dispensing liquid the dispenser should be cleaned.

IX. Repair Service

1. All OPTIFIX® dispensers are made from highest quality materials and are almost maintenance-free and reliable. If a repair should be necessary, it must be done by the manufacturer or an authorized service station. This does not refer to parts available as spare parts (finger guard with piston, volume setting system, discharge tube set, tubes, adaptors).
2. Before returning the dispenser for repair please make sure that it is perfectly cleaned to avoid any contamination and include a certificate of decontamination.
3. Please use the dispenser following to the operating manual to be entitled to the warranty. Repairs must be carried out by the manufacturer or authorized service stations.

X. Technical Data

1. Volume Range - Subdivisions - Accuracy

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*not acc. to ISO 8655/5

The OPTIFIX® dispensers are calibrated to deliver (EX).

The data for accuracy and precision are established according to the German calibration ordinance (PTB) at a room temperature of +20° C (+68° F) by using distilled water. When using chemicals with a different viscosity or density recalibration on an analytical balance is needed.

2. Packing

OPTIFIX® dispenser 2ml, 5ml, 10ml and 30 ml:
single carton, 25,5 x 15,3 x 7,0 cm.

OPTIFIX® dispenser 50ml and 100ml: single carton, 29 x 16,5 x 8 cm.

OPTIFIX® dispenser 200ml, 300ml and 500ml: single carton with foam insert,
46,2 x 18,5 x 14,5 cm.

3. Weight, Dimensions

Volume range in ml	Weight in g appr.	Height in mm
0,4 - 2,0	225	150
1,0 - 5,0	240	165
2,0 - 10,0	300	180
6,0 - 30,0	485	235
10,0 - 50,0	605	245
20,0 - 100,0	800	275
40,0 - 200,0	1500	365
60,0 - 300,0	1800	380
100,0 - 500,0	1900	435

XI. Equipment


All models come in a break-proof single carton with 1 intake tube and 1 discharge tube set. The models BASIC and SOLVENT are additionally equipped with 2 PP adapters, the models SAFETY, SAFETY S and HF are equipped with 2 PTFE adapters.

XII. Warranty


The OPTIFIX® dispenser is probably the most reliable bottle top dispenser in the world. The manufacturer guarantees manufacture and original parts for a 12 month period but excluding all claims for damages, conversion and other claims. Under this warranty all parts that are defective due to faulty materials or production will be replaced. The manufacturer is not responsible for any damages or problems caused by transport or improper handling. All claims must be reported to the supplier. After a warranty repair the warranty period will not be extended. This warranty is not transferrable and void if repairs are carried out by unauthorized service stations.

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All rights reserved, especially the rights to illustrations and original texts in whole or in part, for any and all reproduction.
Copyright 09/2009, by Poulten & Graf GmbH, D-97877 Wertheim